

Delphic AP System Requirements

	Server Requirements	
Hardware	Intel X86-64 compatible	
Processor	Minimum: 2 x CPU	Recommended: 4 x CPU
Memory (RAM)	Minimum: 8GB	Recommended: 16GB
Disk	200GB	
Operating System	Windows Server 2016 Standard Edition	
Server Application	IBM Domino Server 9.0 32 bit	

¹ A virtual server under VMware is supported.

	Client Hardware and Software Requirements	
Hardware	Intel X86-64 compatible	
Processor	Minimum: 1 x CPU	Recommended: 2 x CPU
Memory (RAM)	Minimum: 4GB	Recommended: 16GB
Disk	Minimum: 40GB	Recommended: 100GB
Display	Core Minimum: 1280 x 1024 or SPW: 1920 x 1080 touch screen	
Operating System	Windows 7 or Windows 10	
Client Application	Core: IBM Notes client 9.0 or SPW: Internet Explorer 11 or higher	